


Web Ambassador Quick Start


Release 10 ■ For referring physicians


Contents


- Accessing patient exams via the Internet 1
- The Web Ambassador workspace 1
- Receiving exams automatically and manually 2
- Marking exams as reviewed 3
- Viewing images and listening to audio 3
- Keyboard shortcuts 4
- Getting Help and Information 5


Accessing patient exams via the Internet

- 
1 **Start Web Ambassador.** Double-click the **Web Ambassador** icon.

- 
2 **Select the facility.** If your facility is not displayed in the **Database** field, click the arrow to select it from the **Database** list.


- 
3 **Log on.** Enter your user name and password.

- 
4 **Click the Exam Server tab.**
 If you need to see exams for patients who are in your care, but you did not refer them to the imaging facility, select the **Search All** checkbox.

- 
5 **Search for the patient.** In the Search box, type the first letters of the patient's last name. View the images and reports.

- 6** **Getting help.** If you have questions or need help, do one of the following:
 - Click the **Getting Started** menu.
 - Press **F1**, or On the **Help** menu, click **Help Topics**.


The Web Ambassador workspace



Labels for the Exam Server tab screenshot:

- Menu Bar
- Toolbar
- Database
- Exam Server tab
- Exams list
- Exams you have viewed are highlighted
- Title Bar
- Search field
- Table headings (click to sort)

The Exam Server tab




Labels for the exam viewing window screenshot:

- Menu Bar
- Toolbar
- Viewing tools
- Image viewer
- Title Bar - with patient name and exam date

The exam viewing window

Receiving exams automatically and manually

Receiving exams automatically



Exam Receive Options dialog box

Why use automatic receive?

By receiving exams automatically, you can reduce the time waiting for an exam to download over an Internet connection.

What can you specify?

When you automatically receive exams, you can specify:

- The exam components you want to receive.
- The types of exams you want to view.
- The type of notification, if any – Play a sound or run a program.

How do you display this dialog box?

- 1 On the Receive menu, click Receive Options.
- 2 Complete the fields in the Exam Receive Options dialog box.

Receiving exams manually



Receive Exams Now dialog box

Why request “receive exams now”?

This allows you to check if any exams are currently available that match the criteria you set for receiving automatically.

Displaying this dialog box

- 1 On the Receive menu, click Receive Exams Now.
- 2 Select one or multiple exams from the list, and click Receive.

Disable Standby or Hibernation mode


Requirements to receive exams automatically:

- Your computer must be powered on.
- Web Ambassador must be open.
- The computer must not be in standby or hibernation mode.

To disable hibernation and standby mode:

- 1 Click Start, point to Settings, and click Control Panel.
- 2 Double-click Power.
- 3 Change System Standby and System Hibernates to Never.

Marking exams as reviewed



Why mark an exam as reviewed?

After you review an exam, you can mark it as reviewed. At a glance, you can see which exams you still need to review.

How do you mark the exam?

- 1 Select an exam in the **ExamServer** tab.
- 2 Right-click the exam and click **Mark as Reviewed** in the shortcut menu.

How do you know if an exam has been reviewed?

- 1 Look in the **Reviewed** column.
- 2 Check if the exam has **Yes** in the column. If **Yes** is displayed, the exam has **Yes** in the column, you have viewed the exam.

Viewing images and listening to audio

Main toolbar In the Image Viewer window

Click:	To take this action:
	Return to the Exam Server tab. A message asks if you want to mark the exam as completed.
	Display Properties.
	Compare exam with prior exams.
	Print/fax a report.
	Listen to a recorded audio summary.
	View a report.
	CAD Open the CAD Structured Reports list to view Structured Reports for the exam.
	Copy the exam to a removable media or network file location.
	Display scanned documents.
	Logout from Web Ambassador.
	Return to the Exam Server tab. Suspends viewing.

Selecting and comparing images

In the Image Viewer

To:	Take this action:
Scroll through images:	Click and hold the left mouse button, and move the mouse up or down. OR Press the Up arrow or Down arrow.
Select a series:	Press F11 .
Display another series:	Press Page Up or Page Down .
Specify the number of images displayed:	


Press:	/	1	2	3	4	6	8	9	0	-	=
# of images per screen:	Form at to fit – (all series)	1	2	3	4	6	8	9	10	20	24

Compare exams:
 Click the **Compare** button to display one or more previous exams for the same patient.

Viewing tools

Toolbox

Click to open the toolbox of toolbar buttons and exam information.




Using shortcut menus

- Click the tool button and position the pointer on the image.
- Right-click. From the shortcut menu, select an option.


























Example









Shortcut menu for Page

Through Series 



Viewing tools buttons

Button	Action
	Full Screen Turns off the menus and standard toolbar so the image displays in the full viewing area.
	Montage Displays or hides the montage.
	Set Image-per-Square A viewing mode that displays one image of a series in each viewing window.
	Set Series-per-Square A viewing mode that displays one series in each viewing window. One image in each series is visible at a time.
	Image Panel Properties Changes the display layout, image layout method, presentation rules, and cursor style.
	Previous Panel <i>In Image per Square mode only.</i> Displays the previous panel for the exam being viewed.
	Next Panel <i>In Image per Square mode only.</i> Displays the next panel for the exam being viewed.
	Previous Series <i>In Image per Square mode only.</i> Displays the previous series for the exam being viewed.
	Next Series <i>In Image per Square mode only.</i> Displays the next series for the exam being viewed.
	Select Series Allows you to select a series to display, rather than paging through an entire exam.
	Select and Drop Swaps the position of one series with another selected series.
	Scroll <i>In Image per Square mode only.</i> Scrolls the images on the screen.
	Page Series <i>In Series per Square mode only.</i> Pages through a stack of images in a series.
	Presentation States Displays Presentation States.
	Text Annotation Writes a text note on the image.
	Number Places a number on the image in a specific location.
	Arrow Annotation Places an arrow on the image.
	Measure Annotation Displays measurements with a line or an angle. Also measures stenosis.
	Window/Level Adjusts the contrast of the image.
	Magnify Magnifies the entire image.
	Mag Box Magnifies a portion of the image. Default setting in the Image per Square display mode.
	Invert Shows the negative of the image.
	Rotate/Flip Mirrors the image or rotates the image.
	Crop Saves the selected portion of the image, and removes the rest of the image.
	Calibrate Image Displays pixels per millimeter calculated from a line drawn over the image. Also sizes catheters.

Button	Action
	ROI Marks a region of interest by inserting a pre-selected shape or by drawing a rectangle.
	Shutter Uses a rectangular or circular mask to temporarily block pixels outside the selected area. Used to shut out unwanted light.
	Cardiac Measurement Annotation Displays cardiac measurements.
	Apply Image Processing Applies an edge enhancement algorithm.
	Roam Image Displays an area within the image, when the image is larger than the viewing window. You can move this viewing area across the image.
	Cine <i>In Series per Square mode only.</i> Shows the images in a continuous loop.
	3D Cursor Targets an anatomical area in the primary series, and in all displayed series.
	Cross-ref Scout Creates a scout inset in a series.

Keyboard shortcuts

From the Help menu, click Keyboard Shortcuts.



Alpha keys

* Press the key again to return to the previous state

Press	To
A	Display annotation or DICOM header information.
B	Display the barcode dialog box.
C	Display in Series-per-monitor mode.
D	Display in Exam Compare mode.
E	Expand the image to fill its window.
F	Flip between the primary and other image. Release to return to the primary image.
G	<i>In Ultrasound Stress Echocardiogram only:</i> Stress Default
H	<i>In Ultrasound Stress Echocardiogram only:</i> Stress Images by Stage.
I	Display the scout insets.
K	Display CAD markers.
L	Loop images in continuous display.
O	Show slice line on orthogonal images.
P	Play a pre-recorded voice clip.
Alt+P	Display the Image Panel Properties
Q	<i>In Series-per-square mode:</i> Switch between the "page through series" tool and the last tool used.

S	Switches between Shrink On and Shrink Off. <ul style="list-style-type: none"> Shrink Off magnifies the image to its full acquired resolution. If the image is larger than the viewing square, the Roam tool automatically opens so you can scroll through the image, moving it up and down, left and right, until you view the entire image. Shrink On reduces the image to fit into the image square. The active tool will be left unchanged unless it is the Roam tool, in which case the tool you were using before selecting Shrink Off will be selected.
T	Display the Tools popup menu.
U	Display User annotations.
W	Switch between Wheel and Standard modes.
X	Display in Image-per-square mode.
Z	Display in Series-per-square mode.

Function keys

Press	To
F1	Display Online Help.
F9	Display primary or comparison exam report.
F10	Copy case or single image to removable media.
Alt+F10	Copy a single image to removable media.
F11	Select another series in the exam.
F12	Display the Edit Series dialog box.

Other keys

Press	To
{	Run cine on all images.
}	Stop running cine.
. (period)	Apply to all images. Caution: This action may affect images that are not currently visible on the display.
End	Display last image.
Esc	Exit exam viewing (the same as clicking Done).
Home	Display first image.
PgDn	In Image-per-square mode: Display the next image in exam. In Series-per-square mode: Scroll through the next series.
Alt + PgDn	In Image-per-square mode: Go to the next series.
Ctrl + PgDn	In Image-per-square mode: Automatically scroll the images in all squares.
Ctrl + PgUp	Stops scrolling.
PgUp	In Image-per-square mode: Display the previous image in exam. In Series-per-square mode: Scroll through the previous series.
Alt + PgUp	In Image-per-square mode: Go to the previous series.
Arrow keys	Move the region slowly
Ctrl + arrow	Move the region fast.
Shift + arrow	Resize the region slowly.
Ctrl + Shift + arrow	Resize the region fast.

Getting Help and Information

<p>Displaying Online Help</p> <p>Display a comprehensive Help program:</p> <ul style="list-style-type: none"> Press F1 or click Help Topics from the Help menu. 	<p>Accessing the Resource Center</p> <p>View and download white papers, user guides, TechNotes, technical specifications and other documents:</p> <ul style="list-style-type: none"> Go to www.drsys.com and click the Resource Center button.
<p>Contacting DR Systems</p> <p>10140 Mesa Rim Road San Diego, CA 92121 http://www.drsys.com</p> <p>Phone: 1-800-794-5955 Fax: 1-858-625-3337</p> <p>E-mail: Sales Sales@drsys.com Info Info@drsys.com Support drsupport@drsys.com</p> <p>Customer Support Online:</p> <ul style="list-style-type: none"> Go to www.drsys.com and click Customer Support Login. 	<p>Copyright and Notice</p> <p>Copyright © 1997-2014 DR Systems, Inc. and its licensors. All rights reserved.</p> <p>Protected by U.S. Patents: http://www.drsys.com/legal/</p>

Template: TCP-000180-H